

Sterisets®

Medical Products

Sterisets® Single Use Steel Instruments



Sterisets International BV
Menhirweg 16 | 5342 LS Oss | The Netherlands
P +31 (0)412 667 755 | E info@sterisets.com

www.sterisets.com

PART OF STI GROUP





Sterisets International BV was established in 1996, as a result of a spinoff from a medical distributing company in The Netherlands. Instead of only buying and selling, the goal was to start our own manufacturing facility in the medical business. We have been growing our activities worldwide, with currently 4 facilities in The Netherlands and in Portugal, which together form the STI group.

Together, all of us at Sterisets International Group (STI) are dedicated to provide you with a wide range of medical products to increase efficiency in the medical field.



Sterisets®

Oss, The Netherlands

This is the headquarter of the Sterisets International Group, with the management, international sales division and logistics department.



Sterisets®

Penafiel, Portugal - Logistics & Offices

With this facility, opened in 2020 across the road from our production facility, we have expanded our warehouse capacity with 2.600 m², and 700 m² office space.



Sterisets®

Penafiel, Portugal - Production & Quality

Our 3.500m² facility, gives space to our cleanroom, injection molding production, warehouse and laboratory. In our 1.000 m², ISO 7 cleanroom, thousands of medical kits and trays are produced daily, under close supervision of our highly qualified quality department. The Prefilled syringes are made in a Cleanroom ISO-5, within the same facility.



Sterisets®

Penafiel, Portugal - Non-sterile Production and Warehouse

At our latest Portuguese facility, the production for the lamination and folding of our drapes and draw sheets and storage of our non-sterile products takes place.



Mora, Portugal - Ethylene Oxid (EO) Sterilization and Laboratory Services

With more than 15 years of experience, we support manufacturers across Europe in ensuring the safety, compliance, and performance of their products through validated sterilization processes and high-quality analytical testing.



MEDIRM®

Mora, Portugal - Production and Packaging

MEDIRM Lda is in the market since 2006, is a Portuguese producer whose commercial activity consists in the production and packaging of medical devices with its own brand or client's brand, acting in several European countries.



Always want to have access to the latest version of our brochures? All brochures and leaflets are available digitally on our website www.sterisets.com. In case you prefer a physical copy of our artwork, Sterisets International Group uses Blue Label paper, made out of 100% recycled paper with FSC quality mark. In addition, the paper is produced without usage of any type of chlorine (ECF).

© 2026 Sterisets International B.V. All rights reserved. All images, graphics, text, and other content contained in this brochure are the property of Sterisets International B.V. and are protected by copyright laws. No part of this content may be reproduced, distributed, or transmitted in any form or by any means without the prior written permission of Sterisets International B.V. Unauthorized use is strictly prohibited and may result in legal action. While every effort has been made to ensure the accuracy of the information contained within this brochure, Sterisets International B.V. reserves the right to correct any mistakes or typing errors and shall not be held responsible for any inaccuracies or omissions. Prices, specifications, packaging and availability are subject to change without notice.

Sterisets®

Medical Products

Sterisets® Single use Instruments

Advantages

Quality: High quality recyclable stainless steel.

Hygienic: Every procedure a new clean instrument. Minimize infections.

Logistics: No more logistic problems with cleaning and re-sterilizing instruments.

Cost: Cost effective. Save time and money.

Efficient: Can be inserted in sets & procedure trays.

Colour coded: No confusion: the blue colour together with the single use logo indicates it is a single use instrument.

► Sterilization method: R EO

Colour Coded

Our colour coded single use instruments are made of AISI-410 recyclable stainless steel.

To accentuate the single use of the instruments, next to the single use logo, we mark the instruments with a 100% non toxic blue paint. Extensive tests have proven the safe use of this paint, test rapport is available.

Our instruments are CE marked.

Limitless Possibilities

A wide range of instruments makes it possible to find the right product for every type of procedure, both for use in our custom made sets and trays as single packed sterile in a convenient dispenser box.

Areas of use

Emergency room: Single packed; Suture sets.

Intensive care: Single packed; Emergency sets like IV line pack and tracheotomy.

General ward: Single packed; Suture removal sets.

Private hospitals: Single packed; Full OR sets including drapes and gowns.

Army: Emergency stock.

Single packed sterile

All instruments can be ordered single packed sterile, in case of sharp instruments, we use a protection card to avoid the package from being damaged. We describe this as "".





Table of contents

Introduction	2
Sterisets® Single use Instruments	3
1. Cannulas	5
2. Curettes	6
3. Dental	8
4. Forceps	10
5. Needle holders	19
Sets with Instruments	22
6. Probes	24
7. Retractors	27
8. Scissors	30
9. Speculums	37
10. Miscellaneous	39



Cannula Button
needle Ø 1,2 mm
button Ø 1,8 mm
length 11,5 cm

REF

2010030
* 2010034



Cannula Button
straight needle
length 10 cm

REF

2010082 Ø 1,2 mm
button Ø 2,5 mm
* 2010084

2010086 Ø 2 mm
button Ø 3,25 mm
* 2010088



Cannula Button,
straight needle
length 8 cm

REF


2010014 Ø 1 mm
* 2010015
2010018 Ø 1,2 mm
* 2010052
2010020 Ø 1,5 mm
* 2010054
2010024 Ø 2 mm
* 2010066
2010022 curved Ø 1,5 mm
* 2010062



Cannula Suction
pool and straight
22,5 cm

REF

2010044
* 2010134



Cannula Suction
Frazier
11 cm

REF

2010074 tip size 7
2010076 tip size 8
2010078 tip size 10
2010080 tip size 12

2. Curettes


* **STERILE REF**



Curette Bone
double ended
7/6 mm
length 21 cm

REF

2020014
* 2020016



Curette bone Barth
double ended
18 cm

REF

2020002
* 2020004



Curette Ear Gross
double ended
blunt hook 5mm
spoon 4 x 7 mm
length 23 cm

REF

2020006
* 2020008



Curette Bone
Volkmann
length 17 cm

REF

2020080 Ø 4,4 mm
* 2020220
2020082 Ø 2,8 mm
* 2020214



Curette Evisceration Bunge small tip \varnothing 8 mm, length 13 cm

REF
2020084
*2020086



Curette Bone Volkmann type A double ended 4 x 6 / 6 x 8 mm length 21,5 cm

REF
2020066
*2020068



Curette Bone Volkmann type C double ended 6 x 8,5 / 8 x 12,5 mm length 21,5 cm

REF
2020010
* 2020012



Curette bone Volkmann type B double ended 5 x 6 / 7 x 10 mm length 21,5 cm

REF
2020122
*2020070



Curette Bone Lupus Besnier 5 x 2 mm length 14,5 cm

REF
2020064
*2020208



Curette Bone length 13 cm

REF
2020074 oval tip size \varnothing 2 mm
2020072 round tip size \varnothing 8 mm



Explorer Dental rigid
double ended
1,75 / 1,5 mm
length 17 cm

REF

2020088
* 2020090



Excavator Dental
double ended
Ø 1,3 mm
length 17 cm

REF


2020076
* 2020078



Dental Filling
remover
double ended
length 16 cm

REF

2030008
* 2030010



Dental Mirror no. 5
Ø 2,4 mm handle
length 17 cm

REF

2030002
* 2030004



Probe dental
single ended
14 cm

REF

2060026
* 2060028



Probe Dental
3 mm
graduated
angled 92°
length 14 cm

REF


2060032
* 2060380



Raspatory/Elevators
straight 15 cm

REF


2020032 Fig. 5S
2020038 Fig. 4S
* 202264
2020034 Fig. 3S
* 2020036



Probe Podotherapy
single ended
Ø 1,4 mm
length 13,5 cm

REF

2060090
* 2060092



Probe Periodontal
pocket Williams
graduated
length 14 cm

REF


2060150
* 2060152



Retractor Obwegeser
curved downwards
42 x 11 mm
length 21,5 cm

REF

2070062
* 2070064




Cutting Walls Back
cutting instrument
length 16,5 cm

REF

2090794
* 2090796

4. Forceps

* **STERILE REF**



Forceps Adson 1 x 2 teeth
length 12 cm

REF
2040046
* 2040690



Forceps Adson
dissecting
length 12 cm

REF
2040042
* 2040688



Forceps Adson micro
1 x 2 teeth
length 12 cm

REF
2040062
* 2040696



Forceps dissecting
Adson micro
with plateau
length 12 cm

REF
2040158
* 2040708



Forceps tissue Adson micro
1 x 2 teeth
length 15 cm

REF
2040234
* 2040236



Forceps Allis intestinal grasping
length 15,5 cm

REF

2040420 2 x 3 teeth
* 2041020
2040460 4 x 5 teeth
* 2040462
2040498 5 x 6 teeth
* 2041022



Forceps Baby Allis
4 x 5 teeth
length 15,5 cm

REF

2040178
* 2041024



Forceps artery Baby Mosquito curved
length 10 cm

REF

2040004
* 2040678
2040002 1 x 2 teeth
* 2040954



Forceps Bengolea artery curved
length 24,5 cm

REF

2040394
* 2040978



Forceps Bozemann uterine dressing S-shape
length 25,5 cm

REF

2040492
* 2041116



Forceps Bozemann uterine dressing curved
long. groove
length 25,5 cm

REF

2040490
* 2041122



Forceps Ophthalmic Colibri with tying platform
length 7,5 cm

REF

2040424



Forceps Ophthalmic Colibri 1 x 2 teeth
tip size 0,12 mm
length 7,5 cm

REF


2040454
* 2040944



Forceps Crile artery curved
length 14 cm

REF

2040150
* 2040152



Forceps Dressing Cheron
serrated
handles angled
with ratchet
length 25 cm

REF

2040212
* 2040214



Forceps Dressing uterine
Collin-Pozzi
length 25 cm

REF

2040274
* 2040276



Forceps Dissecting
length 18 cm

REF

2040228
* 2040836
2040676 delicate 15,5 cm



Forceps Dissecting

REF

2040342 11,5 cm
* 2040920
2040034 13 cm
* 2040036
2040058 14,5 cm
* 2040694
2040536 16 cm



Forceps Allis 5 x 6 teeth
intestinal grasping
length 25,5 cm

REF


2040146
* 2040706



Forceps Dressing
Wilde serrated
straight
knee bent
length 11,5 cm

REF

2040516
* 2040518



Forceps Dissecting medium
serrated



REF

2040130 11,5 cm
* 2040870
2040312 25 cm
* 2040314




Forceps Ear Wullsten
length 8 cm

REF
2040526
* 2040936



Forceps Gillies tissue
1 x 2 teeth
length 15,5 cm

REF
2040054
* 2040056



Forceps Dissecting Gillies
length 15,5 cm

REF
2040316
* 2040654




Forceps Dressing Gross
serrated
curved
length 20 cm

REF
2040384
* 2040386



Forceps Dressing Maier
serrated
straight
length 25 cm

REF
2040122
* 2040124



Forceps Hemostatic micro
Mosquito Halsted
serrated
straight
length 12,5 cm

REF
2040184
* 2040186



Forceps artery Micro Halsted
straight 1 x 2 teeth
length 12,5 cm

REF
2040188
* 2040190



Forceps Halsted Mosquito
curved
length 12,5 cm

REF
2040012
* 2040682





Forceps Hemostatic Mosquito
Halsted delicate
straight
length 12,5 cm

REF
2040016
* 2040018



Forceps Tissue, delicate,
straight, 1 x 2 teeth
length 10 cm

REF
2040366
* 2040590



Forceps Dressing Iris
delicate
serrated
curved
length 10 cm

REF
2040428
* 2040430



Forceps Hemostatic Kocher
delicate
serrated,
straight 1 x 2 teeth
length 14 cm

REF
2040084
* 2040700



Forceps Hemostatic
Kocher-Ochsner
serrated
straight 1 x 2 teeth
length 18,5 cm

REF
2040326
* 2040328



Forceps Catheter Introducing
Magill for infants
knee bent

REF
2040162 15 cm
* 2040164
2040166 17 cm
* 2040168
2040170 20 cm for children
* 2040172
2040174 25 cm for adults
* 2040710



Forceps Hemostatic baby
Mosquito Hartmann
fine curved 1 x 2 teeth
length 10 cm

REF
2040020
* 2040022



Forceps Dressing Atraumatic
DeBakey
straight
jaw 2 mm
length 20 cm

REF
2040144
* 2040850



Forceps Debaquey tip 2 mm
straight
length 15,5 cm

REF
2040024
* 2040026



Forceps dissecting Mc Indoe
length 15 cm

REF
2040050
* 2040692



Forceps Micro JFC-7
curved
length 12 cm

REF
2040264
* 2041040



Forceps tissue 1 x 2 teeth
standard
length 14,5 cm

REF
2040066
* 2040698



Forceps tissue with teeth
length 14,5 cm

REF
2040038
* 2040686



Forceps sponge
Rampley serrated

REF
2040100 18 cm
* 2040102
2040104 25 cm
* 2040106



Forceps tissue
medium 1 x 2 teeth
length 18 cm

REF
2040402
* 2040882




Forceps clamp tubing
length 14,5 cm

REF
2040080
* 2040082



Forceps clamp tubing
smooth jaw
length 16,5 cm

REF
2040028
* 2040652



Forceps tissue delicate
1 x 2 teeth
length 11 cm

REF
2040116



Forceps clamp tubing
length 20 cm

REF
2040338
* 2040340



Forceps splinter
Feilchenfeld

REF
2040092 9,5 cm
* 2040094
2040096 11 cm
* 2040704



Forceps sponge Maier
curved ratchet
length 26 cm

REF
2040380
* 2040382



Forceps dressing straight
with ratchet
length 14,5 cm

REF
2040132
* 2041080



Forceps Stieglitz
sharp point
length 13,5 cm

REF
2040118
* 2040120



Forceps Semken dissecting
delicate
length 12,5 cm

REF
2040230
* 2040908



Forceps Dressing delicate
serrated
straight
length 10,5 cm

REF

2040134
* 2040136



Forceps Hemostatic Pean
delicate
serrated
length 14 cm

REF

2040076
* 2040078



Forceps Dressing medium
serrated
length 14,5 cm

REF

2040138
* 2040646



Forceps Rochester Pean
curved
length 18,5 cm

REF

2040154
* 2040156



Forceps artery Rochester Pean
curved and serrated

REF

2040286 20 cm
* 2040296
2040192 22,5 cm
* 2040194
2040294 26 cm
* 2040712



Forceps Foreign Body
Removal Pelkmann
straight
screw lock
length 26,5 cm

REF


2040504
* 2040506



Forceps artery
Rochester Pean straight

REF


2040140 16 cm
* 2040142
2040070 18,5 cm
* 2040072



Forceps uterine tenaculum
Schröder
length 25 cm

REF


2040204
* 2040206



Forceps Schumacher Biopsy
working length
length 24 cm

REF


2040494
* 2040650



Forceps Spencer Wells
artery straight
length 20 cm

REF


2040352
* 2040354



Forceps uterine Pozzi
fine pattern
length 25 cm

REF

2040404
* 2041128



Forceps Dressing Troeltsch,
serrated
straight
knee bent
length 12 cm

REF

2040126
* 2040128



Forceps Tracheal Dilator
Trousseau
length 14 cm

REF

2040108
* 2040110





Needle holder Baby
Crile-Wood
length 15 cm

REF
2050016
* 2050208



Needle holder Crile-Murray
serrated
length 15 cm

REF
2050144
* 2050146




Needle holder Crile-Wood
serrated
length 18 cm

REF
2050124
* 2050272



Needle holder Debaquey
serrated
length 15 cm

REF
2050114
* 2050116





Needle holder Baumgartner
length 14,5 cm

REF
2050010
* 2050012



Needle holder Halsey
length 13 cm

REF
2050074
* 2050294



Needle holder Crile-Wood
length 14,5 cm

REF
2050048
* 2050049



Needle holder Lawrence
length 15 cm

REF
2050120
* 2050186



Needle holder Crile-Wood
fine pattern, serrated
length 15 cm

REF
2050014
* 2050282



Needle holder Wertheim
serrated, straight
bayonet shaped
length 20 cm

REF
2050080
* 2050364



Needle holder Derf
length 12,5 cm

REF
2050122
* 2050290



Needle holder Ryder
smooth jaw
length 14 cm

REF
2050042
* 2050340



Needle holder Hegar
length 13 cm

REF
2050026
* 2050028



Needle holder
Micro Arret
length 16 cm

REF
2050064
* 2050336



Needle holder Mathieu
length 14 cm

REF
2050034
* 2050298



Needle holder
smooth jaw
length 12,5 cm

REF
2050136
* 2050138



Needle holder Mayo-Hegar
serrated
length 15 cm

REF
2050088
* 2050090



Needle holder
Mayo-Hegar

REF
2050002 15 cm
* 2050206
2050022 16 cm
* 2050024
2050036 18 cm
* 2050038
2050050 20,5 cm
* 2050318



Needle holder Ryder
with pattern
14 cm

REF
2050044
* 2050046

Sterisets®

Medical Products

Single use Instruments

Sets and Trays

Build your customized medical sets and trays containing Sterisets colour coded single use instruments.

The instruments are manufactured by an ISO 9001 & ISO13485 certified company and meet the international standards.


You can verify some examples of our sets that contain our instruments.





ANESTHESIA SETS

EXAMPLE - REF. 1011016



EYE SETS

EXAMPLE - REF. 1060002



NEFROLOGY SETS

EXAMPLE - REF. 1300040



PLEURAL SETS

EXAMPLE - REF. 1300014



SUTURE SETS

EXAMPLE - REF. 1050004




SUTURE REMOVAL SETS

EXAMPLE - REF. 1050020

6. Probes



* **STERILE REF**



Probe Gall Duct malleable
straight
double ended Ø 3,5 mm
length 25 cm

REF



2060002
* 2060388



Probe Director
straight
length 14,5 cm

REF


2060020
* 2060022




Probe Pointed
Ø 1,4 mm
length 14 cm

REF

2060082






Probe Gall Duct malleable
straight
double ended Ø 3,5 mm
length 25 cm

REF

2060002
* 2060388



Probe Myrtle Leaf malleable
straight
Ø 1,5 mm
with tip
length 14,5 cm

REF


2060012
* 2060014



Probe Uterine Sims malleable
graduated
curved
length 33 cm

REF

2060038
* 2060040



Probe Director straight
length 14,5 cm

REF


2060020
* 2060022



Probe Foreign
body Zaufal
length 16 cm


REF

2060112
* 2060114




Probe with eye
length 14,5 cm

REF
2060188



Probe Bowman malleable,
straight
double ended \varnothing 0,7/0,8 mm
length 13 cm

REF
2060008
* 2060010



Probe Wilder
lacrimal dilator
double ended, size
000/0000
length 12,5 cm

REF
2060030
* 2060400





Retractor Wound baby
Senn-Miller
3 prongs 8 x 7 / 22 x 7 mm
length 15,5 cm

REF

2070046
* 2070048



Retractor Small curve
blunt
2 prongs 4,5 x 4,2 mm
length 16,5 cm

REF

2070058
* 2070060



Retractor Brunner

REF

2070092 8 x 2 x 27,5 cm
* 2070262
2070094 10 x 2,5 x 27,5 cm
* 2070256




Retractor Lid Desmarres


REF

2070002 depth: 9 x 9 mm
length: 14 cm
* 2070004

2070006 depth: 9 x 12 mm
length: 16 cm
* 2070008

2070010 depth: 11 x 14 mm
length: 16 cm
* 2070012






Retractor eye Barraquer
straight edges

REF
2070034 length: 4 cm
* 2070286
2070038 length: 4,5 cm
* 2070040



Retractor Langenbeck
3,3 x 1,4 cm
length 21 cm

REF
2070028
* 2070030



Retractor Wound Logan
modified Reill
self retaining
3 x 3 prongs
length 5 cm

REF
2070080
* 2070322



Retractor Wound Ragnell
double ended
5 x 3 / 14 x 6,3 mm
length 15 cm

REF
2070014
* 2070218




Retractor Weitlaner
self retaining
sharp/deep
3x4 prongs
length 14 cm

REF
2070170
* 2070218



Retractor Kocher-Langenbeck
single ended
41 x 11 mm
length 21,5 cm

REF
2070114
* 2070228



Retractor Skin Gillies
large
sharp
1 prong
length 18 cm

REF
2070032
* 2070348



Retractor Weitlaner
self retaining
sharp
3 x 4 prongs
length 13 cm


REF
2070104
* 2070106



Retractor Weitlaner
self retaining blunt
length 11 cm

REF

2070018
* 2070220



Retractor Wound / Trachea
sharp large curve
1 prong
length 16,5 cm

REF

2070054
* 2070056



Retractor Mannerfelt
small curve
sharp
2 prongs 5 x 4 mm
length 15,5 cm

REF

2070076



Retractor Volkmann
4 prongs
depth 9 x 19 mm
length 22 cm

REF

2070024
* 2070026



8. Scissors

* **STERILE REF**



Scissors Braun
Stadler
episiotomy
length 14,5 cm

REF

2080086
* 2080088



Scissors curved
sharp/blunt
length 13 cm

REF

2080064
* 2080066



Scissors suture
sharp/blunt
length 13 cm

REF

2080060
* 2080062



Scissors suture
Spencer
length 11,5 cm

REF

2080334
* 2080336



Scissors with
probe at end
length 14 cm

REF

2080090
* 2080092



Scissors
bandage
Lister
length 14 cm

REF

2080014
* 2080448



Scissors bandage
Lister
with 1 large ring
length 18,5 cm

REF


2080284
* 2080286



Scissors curved
sharp/sharp
length 14,5 cm

REF

2080046
* 2080048



Scissors
dissecting Mayo
Stille straight
length 17 cm

REF

2080208
* 2080210



Scissors
dissecting Mayo
curved

REF

2080054 15,5 cm
* 2080512
2080427 17 cm beveled blade



Scissors
dissecting Mayo straight

REF

2080206 15 cm chamfered
* 2080524
2080272 23 cm
* 2080274



Scissors
suture
straight
blunt/sharp
length 18,5 cm

REF

2080146
* 2080460



Scissors
bandage
Lister
serated
length 16 cm

REF
2080082
* 2080084



Scissors
dissecting
curved
Metzenbaum
Baby
length 14,5 cm

REF
2080138
* 2080140



Scissors
dissecting Kilner
length 15 cm

REF
2080002
* 2080004



Scissors
dissecting
Stevens
curved
length 11 cm

REF
2080430



Scissors
dissecting
Mayo
blunt/blunt
length 14 cm

REF
2080306



Scissors Iris
delicate
curved
length 11 cm

REF
2080108
* 2080626



Scissors
dissecting
curved
Metzenbaum
length 14,5 cm

REF
2080050
* 2080052



Scissors Iris
curved
length 11 cm

REF
2080068
* 2080456



Scissors
dissecting
curved
Metzenbaum
length 18 cm

REF

2080094
* 2080458



Scissors
dissecting
curved
Metzenbaum
Durotip
length 18 cm

REF

2080178
* 2080572



Scissors dissecting
Tenotomy Stevens
curved
length 10 cm

REF

2080200
* 2080202



Scissors Iris
delicate
curved
serrated
length 11 cm

REF

2080098
* 2080100



Scissors Iris
and ligature
curved
length 11,5 cm

REF

2080128
* 2080130



Scissors Iris
and ligature
sharp/sharp
straight
length 11,5 cm

REF

2080124
* 2080126



Scissors Iris curved
sharp/sharp

REF

2080122 11 cm supercut
* 2080642
2080106 13 cm
* 2080708



Scissors Iris
delicate
sharp/sharp
straight
length 11 cm

REF

2080006
* 2040684



Scissors Iris
sharp/blunt
straight
length 11 cm

REF

2080308
* 2080636



Scissors Iris
sharp/sharp
straight
length 10,5 cm

REF

2080218



Scissors suture
Littauer
length 13,5 cm

REF

2080078
* 2080080



Scissors baby
Metzenbaum
delicate
curved

REF

2080034 14,5 cm
* 2080454
2080444 11,5 cm



Scissors baby
Metzenbaum
delicate
straight

REF

2080030 14,5 cm
* 2080452
2080446 18 cm



Scissors
Strabismus
curved
length 11,5 cm

REF

2080074
* 2080076



Scissors nail
Konig
sharp/blunt
length 15 cm

REF

2080244
* 2080246



Scissors suture
straight
sharp/sharp
length 14 cm

REF

2080010
* 2080470



Scissors Dissecting
Baby Metzenbaum
curved
wave cut
length 14,5 cm

REF
2080228
* 2080506



Scissors Suture
Spencer
straight
notch at the tip
length 9 cm

REF
2080042
* 2080044



Scissors straight
sharp/blunt

REF
2080026 14,5 cm
* 2080450
2080270 25 cm



Scissors
Umbilical Cord
blunt/blunt
curved sideways
length 13,5 cm

REF
2080164
* 2080166



Scissors suture
sharp/blunt
curved
length 11,5 cm

REF
2080312



Scissors uterine
Sims
blunt/blunt
curved
length 20 cm

REF
2080038
* 2080040



Scissors Suture
Cooper
straight
blunt/blunt
length 14,5 cm

REF
2080114
* 2080116



Scissors Suture
curved
sharp/blunt
length 14,5 cm

REF
2080056
* 2080058



Scissors Suture
straight
sharp/sharp
length 13 cm

REF

2080132
* 2080134



Scissors Delicate
Dissecting Baby
Metzenbaum
curved
length 14,5 cm

REF

2080300
* 2080140



Scissors Micro
curved
sharp/sharp
length 14 cm

REF

2080340
* 2080660



Scissors Micro
sharp/sharp
straight
flat serrated handle
length 11 cm

REF


2080110
* 2080112



Scissors
Stitch-Westcott
curved
blunt/blunt
length 12 cm

REF

2080236
* 2080748



Scissors Capsulotomy
Vannas
curved
flat serrated handle
length 10 cm

REF


2080102
* 2080104



Scissors Micro
straight
sharp/sharp
flat serrated handle
length 11 cm

REF


2080330
* 2080332



Speculum 2 parts
length 19,5 and
23 cm

REF


2090002
* 2090136



Speculum eye Lieberman
V-shaped wire blades
length 8 cm

REF

2090040
* 2090140



Speculum eye Barraquer
15 mm solid blades
length 4 cm

REF

2090044
* 2090046




Speculum ear
Hartmann

REF

2090094 3,2 mm
* 2090096
2090098 3,5 mm
* 2090137
2090082 4,6 mm
* 2090084
2090100 5 mm
* 2090138
2090102 6,5 mm
* 2090104






Speculum nasal
Thudichum
length 7 cm

REF

- 2090004 blade 8 x 12 mm
- * 2090006
- 2090008 blade 12 x 26 mm
- * 2090010



Speculum nose Killian
length 14 cm

REF

- 2090012 blade 56 x 7 mm
- * 2090014
- 2090016 blade 75 x 7 mm
- * 2090150



Speculum Nasal
Hartmann-Halle
short
blade 25 mm
length 15 cm

REF


- 2090018
- * 2090143



Speculum vaginal Cusco

REF


- 2090020 blade 80 x 24 mm
- * 2090022
- 2090024 blade 80 x 32 mm
- * 2090026
- 2090056 blade 92 x 25 mm
- * 2090160
- 2090072 blade 90 x 36 mm
- * 2090160
- 2090058 blade 100 x 27 mm
- * 2090060
- 2090062 blade 110 x 30 mm
- * 2090159



Speculum vaginal
Kristeller
lower blades

REF

- 2090064 70 x 27 mm - length 18 cm
- * 2090169
- 2090066 90 x 36 mm - length 18,5 cm
- * 2090172





Bone Rongeur Beyer
double action
curved
length 18 cm

REF

2040524



Towel Clamp for paper drape
length 15 cm

REF

2090818

* 2090792



Towel Backhaus
curved

REF

2090782 19 cm

* 2090784

2090786 11 cm

* 2090788

2090790 13,5 cm

* 2090890



Cutter Nail
concave cutting edge
with lock
length 14,5 cm

REF

2040448


* 2040450



Gouge Eye foreign body
tip size 1,5 x 15 mm
length 12 cm

REF


2090734
* 2090872



Hook Skin Gillies small
sharp 1 prong
length 18 cm

REF

2070068
* 2070070



Loop ENT/Ear, Billeau
small, Ø 2 mm
length 16 cm

REF

2060016
* 2060018



Suction tube Style M
tip size 3 mm

REF

2010042
* 2010136



Needle Iowa,
trumper shaped
length 14 cm

REF


2050060
* 2050062



Ear hook Lucae, probe end
angled 90° 1 prong
length 14 cm

REF

2070134
* 2070280



Knife Dura Halle
length 20 cm

REF

2090506
* 2090524



Knife Amputation Liston
blade size 13 cm


REF

2090512



Laryngoscope Miller,
blade no. 0
with plastic connector

REF
2090724




Luer None Rongeur
Blumenthal, curved
length 15,5 cm

REF
2040528
* 2040530



Laryngoscope Macintosh
blade with plastic
connector

REF
2090730 N° 0 8,2 cm
* 2090736
2090738 N° 1 9,2 cm
* 2090740
2090744 N° 2 11,1 cm
* 2090746




Discision Walton
foreign body
blunt Ø 1,4 mm
length 12 cm

REF
2050110
* 2050112



Nerve hook Dandy
blunt tip 8 mm
length 22 cm

REF
2070128
* 2070330



Nerve Adson, blunt
angled 90°, 1 prong
length 19 cm

REF
2070074



Needle Pudental
length 14 cm

REF
2050066
* 2050068



Raspatory Bone
slightly curved, sharp
length 16,5 cm

REF
2090800
* 2090832



Tooth Flag
smooth,
length 15 cm

REF
2040180
* 2041090



Spatula Urethra testing
length 18,5 cm

REF
2060086
* 2060088



Scalpel No. 3, straight
length 12,5 cm

REF
2090502
* 2090524



Staple remover
length 10 cm

REF
2090200
* 2090202



Scratch needle
length 3,7 cm

REF
2090520
* 2090528



Scalpel
length 11,5 cm

REF
2090514
* 2090516





Sterisets®

Medical Products

Sterisets International BV

Menhirweg 16 | 5342 LS Oss | The Netherlands

P +31 (0)412 667 755 | E info@sterisets.com

www.sterisets.com

PART OF STI GROUP

Sti

Sterisets International