

Raw Identification	
Reference	85000
Name	Blanket warming unit
Sterilization	NA
Supplier Identification	
Model	WB1000F

Specifications
Electrical
Power supply: 220~240 V, 50/60 Hz, max. 4.5A
Protection: Over-heating (visual and audio alarm) + automatic shutdown
Power consumption: 500W (1.000W peak)
Technical
Filter type: 0,2µm HEPA filter with efficiency 99,97%
Classification
Device classification: Class IIa
Standards: EN 60601-1-2, EN 60601-1
Performance
Temperature options: Low 32°C / Med 38°C / High 43°C
Temperature accuracy: ±2°C
Temperature sensors: 2
Air Flow
Low fan: 13.21 L/s
High fan: 16.52 L/s
Noise Levels
Low fan: 51.3
High fan: 54.7
Physical
Dimensions: H373 x W223 x D172 mm
Weight: 2,3 kg
Power cord: 4,5m (US, EU & UK plugs are available)
Installation: The unit can be attached to the IV stand or bed rail
Time: LED indicator for replacing filter
Hose: 1,8m length, flexible



Packaging content:

- 1 : Sterisets® Blanket Warming Unit
- 2 : I.V. Pole bracket
- 3 : Power cord 4,5m
- 4 : Power cord safety bracket
- 5 : Air hose securing clip
- 6 : Bed rail bracket
- 7 : Reinforced air hose 1,8m

